



2
Lifting
columns

+

2
Linear
actuators




ALARA Y

Electromechanic Pediatric Bed



Actuators With Min
IPX4 Protection
For Quiet Operation

Ø150 mm 360 ° Swivel
Castors With Central
and Directional Lock



Hand Control Unit
Nurse Panel With
Function Locking Feature



Control Box
Built-in Battery





Manual CPR Arm


Mechanical Angle Indicator
on Siderails (trend angle)



Height Adjustable Iv Pole




ABS Plastic Bottom Chassis




47-77cm Height Movement



12° rendelenburg/Antitrendelnburg



70° Backrest Movement



30° Legrest Movement

Width (Over all)	Lenght (Over all)	Height Min.	Height Max.	Lying Surface	Mattress	Backrest Angle	Legrest Angle	Trendelenburg Angle	Antitrendelenburg Angle	Safety Working Capacity	Operating Voltage	Number of Actuators
103,5 cm	180 cm	47 cm	77 cm	78x162 cm	78x160x12 cm 28DNS	70°	30°	12°	12°	230 kg	220 – 240 V AC	4

4-piece lying surface
made of HPL material

75 mm (± 5 mm)
Plastic protective
bumpers against impact

Nano ceramic coated
electrostatic powder
painted metal parts



Foldable side rails
and removable head/foot boards
made of Polypropylene (PP)

ABS Iv Pole Holder

ABS Plastic
Bottom Chassis
Cover

Bottom chassis profiles:
50x50x3mm
30x50x2mm

Upper chassis profiles:
30x50x2mm
40x80x2mm

Material Features