



Autoclave Machine

(23 Liter)

Usage:

The table top Class B autoclaves are with safe auto-controlled, which are designed for dental clinics, hospitals, laboratories, school labs, scientific research institutions, mostly suit for sterilization of surgical instruments, dental instruments, ophthalmic instruments, glassware medicine, culture media, fabrics, etc.

Main Features:

1. Quality 304 stainless steel inner tank, sealed cover.
2. Microcomputer control.
3. Safety interlocking device.
4. Automatic water feed.
5. With vacuum drying function.
6. With BD test, HELIX test, VACUUM test.
7. LCD display with touch buttons.
8. With standard test interface.
9. European Class B standard.
10. With automatic protection function: over temperature protection, over pressure protection, low water level protection, anti-dry burning, and safety door lock system.
11. 3-times pre-vacuum.
12. Buzzer reminder after sterilization, automatic stop running.
13. Automatic discharge of cold air, automatic steam exhaust after sterilization.
14. Built-in steam generator to provide saturated steam quickly.
15. Standard with Printer, USB interface.

Main Parameters:

Model	HY-TS23B
Volume	23L
Power supply	2KW
Voltage	220~240V, 50/60Hz
Design pressure	0.25Mpa
Design temperature	139°C
Rated working pressure	0.22Mpa
Rated working temperature	134°C
Chamber volume	Ø250*452mm
Overall dimension	480*690*450mm
Packing size	760*560*510mm
G.W./N.W.	60kg/53kg