



Medikal 2000

Technology for Health


3
Linear
actuators



DC-03

Blood Donor Chair

Colour
Options





Actuators With Min
IPX4 Protection
For Quiet Operation

Hand Control Unit



Control Units (Standard)
Battery (Optional)

Control units pictures are representative.



70° Backrest Movement



80° Legrest Movement



67-93cm Height Movement

Width (Over all)	Lenght (Over all)	Height Min.	Height Max.	Backrest Surface	Sitting Surface	Backrest Angle	Legrest Angle	Safety Working Capacity	Operating Voltage	Number of Actuators
95 cm	185 cm	67 cm	93 cm	60x80 cm	60x50 cm	70°	80°	230 kg	220 – 240 V AC	3

