



SonoTrax
Ultrasonic Pocket Doppler

About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- *In-Vitro* Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



Global Headquarters:
Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen
518122 P.R. China | +86.755.26898326 | www.edan.com | info@edan.com
U.S. and Canada inquiries:
EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126
+1.858.750.3066 | www.edandiagnosics.com | edan-info@edandiagnosics.com

© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.



ENG-OB & GYN-SONOTRAX
V3.4-20230206

SonoTrax

Ultrasonic Pocket Doppler



SonoTrax

Ultrasonic Pocket Doppler

SonoTrax Series Ultrasonic Pocket Dopplers' superior sensitivity, high-fidelity sound, recording option and optional rechargeable battery ensure the doppler family to meet the routine examination requirements of obstetricians and midwives.

Specifications

- Display: 45 x 25 mm LCD display
- Heart Rate Range: 50~210 bpm
- Resolution: 1 bpm
- Recording Time: 240 s
- Audio Output Power: 1W
- Auto Power-off: 1 minute with no signal or operation
- Size: 34 (D) x 89 (W) x 141 (H) mm
- Weight: about 300 g (including the battery)



Carry Case



Rolling Stand



SonoTrax II Pro



Compact Battery Charger



Backlit LCD



Built-in Recorder



Interchangeable Waterproof Probe

- Easy-to-use
- Compact design
- Support 1.5 V AA battery
- Built-in audio recorder
- Large backlit LCD display
- Accurate fetal heart detection with high-fidelity sound
- Support 2 MHz, 3 MHz, 4 MHz, 5 MHz, 8 MHz waterproof probes
- Automatic and manual counting of FHR
- Interchangeable waterproof probe
- Automatic power-off



SonoTrax Vascular



SonoTrax Lite



SonoTrax Basic



SonoTrax Basic A



SonoTrax Pro



SonoTrax II



SonoTrax II Pro