



**HYS-01
Patient Wash
Stretcher
(Hydraulic)
STRETCHER**

- Nano ceramic coated electrostatic powder painted metal parts
- 1 sections lying surface made of HPL
- Easy foldable side rails made of metal
- Custom liquid proof pool mattress
- Trendelenburg by gas spring
- 360° turning casters (150 mm) with central and directional lock
- ABS plastic covered lower chassis
- Hydraulically controlled height adjustment
- Bumpers on four corners
- Electrostatic powder coated pushing bars



HYS-01
Patient Wash
(Hydraulic)
STRETCHER



HYS-01
Patient Wash
(Hydraulic)
STRETCHER



**180° Opening Side
Rail System**



**HEIGHT ADJUST.
POSITION**



**TRENDELENBURG
POSITION**



General Features

Width (Over all)	Lenght (Over all)	Height Min	Height Max	Lying Surface			Trendelenburg Angle		Safety Working Capacity
84 cm	208,5 cm	54,5 cm	83 cm	60*193,5 cm			12		230 kgs



HYS-01
Patient Wash
(Hydraulic)
STRETCHER