

IKONIUM 4S

Electromechanic Hospital Bed

4
Linear
actuators



Colour
Options





Actuators With Min
IPX4 Protection
For Quiet Operation

360 ° Sentral Brake
Castors




Hand Control Unit



Control Box



Manual CPR Arm



Mechanical Angle Indicator
on Siderails (backrest- trend angle)



Height Adjustable Iv Pole



Nurse Control Panel



Control units pictures are representative.



Width (Over all)	Lenght (Over all)	Height Min.	Height Max.	Lying Surface	Mattress	Backrest Angle	Legrest Angle	Trendelenburg Angle	Antitrendelenburg Angle	Safety Working Capacity	Operating Voltage	Number of Actuators
109 cm	214 cm	47,5 cm	85 cm	86x192 cm	85x190x12 cm 28DNS	77°	40°	18	18	250 kg	220 – 240 V AC	4

