

1
Lifting
columns

+

2
Linear
actuators



Colour
Options



M 2 - J X
Gynecology Examination Table



Actuators With Min
IPX4 Protection
For Quiet Operation

Control Box
Built-in Battery



Remote Control



Leg Extension



Squatting bar

Waste Basin




Calf Support




Optional Foot Control

Control units pictures are representative.



15° Trendelenburg / 15° Antitrendelenburg



63-105 cm Height Movement



60° Backrest Movement

Width (Over all)	Lenght (Over all)	Height Min.	Height Max.	Lying Surface	Sitting Surface	Backrest Angle	Trendelenburg Angle	Antitrendelenburg Angle	Safety Working Capacity	Operating Voltage	Number of Actuators
66 cm	167 cm	63 cm	105 cm	65x81 cm	52.5x65 cm	60°	15°	15°	230 kg	220 – 240 V AC	3

Two-pieces
upholstery with
polyurethan filling

Anticorrosive
chromium-nickel
covered accessories.

Ability to fixing
on the floor



Bottom chassis profiles:
30x50x2mm

Upper chassis profiles:
20x40x2mm

ABS Plastic
Bottom Chassis
Cover

Nano ceramic coated
electrostatic powder
painted metal parts

Material Features