



Colour
Options



MA - 03
Emergency Trolley / Crash Cart



Disposal bin with lid
made of ABS material



Seal system for
central locking drawers

Telescopic rails for
easy opening and
closing drawers



90° turning defibrillator holder



Nurse writing table

IV Pole with telescopic adjustment



Oxygen tube holder
made of polypropylene



Optional
Key Locking
Feature



Side panels and top table
made of polypropylen material

Main frame made of
durable aluminum material

360 ° rotatable
2 castors with brakes
and 2 non-brake wheels

5 drawers with following dimensions:

2 pcs 445x545x70mm,
2 pcs 445x545x160mm,
1 pc 445x545x250mm

Width	Lenght	Height
79 cm	61 cm	110 cm