

Portable Oxygen Concentrator with battery



OLV-C1 portable oxygen concentrator with Battery
is the most lightest, smallest and quietest oxygen
concentrator for travel and ambulatory oxygen needs.



Zhengzhou Olive Electronic Technology Co.,Ltd.
Tel: 0086 371 8608 1306 Fax: 0086 371 6906 7368

www.zzolive.com
E-mail: info@zzolive.com

FEATURES

1. Continuous flow, 1-5L/min adjustable
2. Light weight: 5kg
3. Three power supplies: 220V/110V, 12V(Car use), Rechargeable Lithium-ion Battery
4. Come with carrying bag and portable trolley. It is more convenient to use at outdoor
5. Low noise, less than 40db(A)
6. Anion function



SPECIFICATIONS

Size	310*165*310mm
Weight	machine 5kg, battery 1kg
Flow Settings	1-5LPM
Purity Oxygen	≥90%±3%(1LPM)
Output Pressure	40Kpa-60Kpa
Sound Level	less than 40db(A)
Power	1. AC Power 220V±22V, 50Hz±1Hz/110V±15V, 60Hz±1Hz 2. DC Power 12V 3. Rechargeable Battery
Power Consumption	90W
Battery Duration	60mins
Battery Recharge Time	4 hours recharge time to achieve 90% capacity