




3
Lifting
columns

+

2
Linear
actuators


PLUS T

ELECTROMECHANIC ICU
HOSPITAL BED WITH THREE COLUMNS



Actuators With Min IPX4
Protection For Quiet
Operation

Accessory Pocket



Hand Control Unit
Nurse Panel
With Function Locking Feature



Control Box
Built-in Battery



Manual CPR Arm



Mechanical Angle Indicator
on Siderails
(backrest-trend angle)



Height Adjustable IV Pole



Ø 150 mm 360 ° Swivel
Central Brake Castors
and ABS Plastic Bottom Chassis



Control units pictures are representative.



File Holder



Oxygen Cylinder Holder




Fifth Caster



Actuators With Min IPX4
Protection For Quiet
Operation



Leg Extension Part



Device Stand


Actuators With Min IPX4
Protection For Quiet
Operation



HPL Lying Surface



18° Trendelenburg/Antitrendelenburg



60-100 cm Height Movement



25° Lateral Tilt Movement



70° Backrest Movement



30° Legrest Movement

Width (Over all)	Lenght (Over all)	Height Min.	Height Max.	Lying Surface	Mattress	Backrest Angle	Legrest Angle	Trendelenburg Angle	Antitrendelenburg Angle	Lateral Tilt Angle	Safety Working Capacity	Operating Voltage	Number of Actuators
103,5 cm	212,5 cm	60 cm	100 cm	85x190 cm	85x190x12 cm 28DNS	70°	30°	18°	18°	25°	230 kg	220 - 240 VAC	5

