



3
Linear
actuators



PLUS A3

Electromechanic Hospital Bed

Colour
Options





Actuators With Min
IPX4 Protection
For Quiet Operation

Ø 125 mm 360 ° swivel
cross lock castors
ABS Plastic Bottom Chassis



Hand Control Unit



Control Box



Manual CPR Arm

Mechanical Angle Indicator
on Siderails (backrest- trend angle)



Height Adjustable Iv Pole




ABS Plastic Bottom Chassis

Control units pictures are representative.



Siderails Control Panel


Underbed Light



Sticker options
in different colours
for Head/Footboards
and side rails



HPL Lying Surface



X-Ray Cassete Holder

Leg Extension Part




Nurse Control




File Holder






43-78,5cm Height Movement



70° Backrest Movement



30° Legrest Movement

Width (Over all)	Lenght (Over all)	Height Min.	Height Max.	Lying Surface	Mattress	Backrest Angle	Legrest Angle	Trendelenburg Angle	Antitrendelenburg Angle	Safety Working Capacity	Operating Voltage	Number of Actuators
103,5 cm	212,5 cm	43 cm	78,5 cm	85x190 cm	85x190x12 cm 28DNS	70°	30°	-	-	230 kg	220 – 240 V AC	3

