


PLUS A7

Electromechanic Hospital Bed

3
Linear
actuators



 **Medikal 2000**
Technology for Health

Colour
Options





Actuators With Min
IPX4 Protection
For Quiet Operation

Ø 125 mm 360 ° swivel
cross lock castors
ABS Plastic Bottom Chassis



Hand Control Unit



Control Box



Manual CPR Arm

Mechanical Angle Indicator
on Siderails (backrest- trend angle)



Height Adjustable Iv Pole




Trendelenburg Arm

Control units pictures are representative.



Siderails Control Panel


Underbed Light



Sticker options
in different colours
for Head/Footboards
and side rails



HPL Lying Surface



X-Ray Cassete Holder


Leg Extension Part




Nurse Control



File Holder




43-84,5cm Height Movement



12° Trendelenburg Movement



70° Backrest Movement



30° Legrest Movement

Width (Over all)	Lenght (Over all)	Height Min.	Height Max.	Lying Surface	Mattress	Backrest Angle	Legrest Angle	Trendelenburg Angle	Antitrendelenburg Angle	Safety Working Capacity	Operating Voltage	Number of Actuators
103,5 cm	212,5 cm	43 cm	84,5 cm	85x190 cm	85x190x12 cm 28DNS	70°	30°	12°	-	230 kg	220 – 240 V AC	3



ABS Iv Pole Holder

75 mm (± 5 mm)
Plastic protective
bumpers against impact

ABS Plastic
Bottom Chassis
Cover

Bottom chassis profiles:
30 x 50 x 2 mm
Rectangular Tube

Upper chassis profiles:
40 x 80 x 2 mm
30 x 50 x 2 mm

Removable and mountable
head and foot boards
made of Polypropylene (PP)

4-piece lying surface
made of ABS material

Nano ceramic coated
electrostatic powder
painted metal parts

Material Features