





3
Linear
actuators



RC - 01
Home Care Bed

Colour
Options





Actuators With Min
IPX4 Protection
For Quiet Operation

Ø150 mm 360 ° Swivel
Central Brake Castors



Hand Conrol Unit



Control Box



Manual CPR Arm


HPL Lying Surface



Height Adjustable Iv Pole



Underbed Light





Two sided seat movements



30° Legrest Movement



70° Backrest Movement



Armrest movement

Width (Over all)	Lenght (Over all)	Height	Lying Surface	Mattress	Backrest Angle	Legrest Angle	Trendelenburg Angle	Antitrendelenburg Angle	Safety Working Capacity	Operating Voltage	Number of Actuators
103,5 cm	214 cm	42 cm	85x190 cm	85x190x12 cm 28DNS	70°	30°	-	-	230 kg	220 – 240 V AC	3

