




T-300
Hydraulic Emergency Stretcher




Height Adjustable Iv Pole

Height Controlling Pedal



Side Rail Lock Arm




ABS Plastic Bottom Chassis
and Central Brake System




Double Shock Absorber
on the Backrest

X-Ray Cassete Holder





90° Backrest Movement



54,5-83 cm Height Movement



12° Trendelenburg

Width	Lenght	Height Min.	Height Max.	Lying Surface	Mattress	Backrest Angle	Trendelenburg Angle	Safety Working Capacity	Trend. Motion Source	Height Motion Source	Backrest Motion Source
84 cm	208,5 cm	54,5 cm	83 cm	60 x 193,5 cm	60x190x8 cm	90°	12°	230 kg	gas powered	hydraulically powered	gas powered

