



**T-700  
HYDRAULIC  
EMERGENCY  
STRETCHER**

- Nano ceramic coated electrostatic powder coated metal parts
- 4-part lying surface made of compact material
- Easy to open and close side rails
- Crash bumpers on all four corners
- Telescopic adjustable IV pole
- Plastic coated foldable push bars
- 360° swivel central locking wheels (200 mm)
- 5.Th Castor (Provides ease of movement in the desired direction on the stretcher)
- Hydraulically controlled height, trendelenburg and reverse trendelenburg adjustment with foot pedals
- Gas spring controlled back and foot adjustment
- Ability to adjust knee movements with gradual mechanism
- Aluminum telescopic columns
- ABS plastic coated sub-chassis
- Oxygen cylinder holder
- Patient restraint belt
- 28 density mattress (66x193x8cm)
- X-Ray Cassette carrier




**T-700  
HYDRAULIC  
EMERGENCY  
STRETCHER**



1

Ayarlanabilir serum askısı  
Height adjustable IV pole



2

Ayak Pedalı ile kolay Kullanım  
Foot Pedal For Easy Use



3

5. Teker  
5.Th Castor



4

ABS Alt Şase ve Merkezi Kilit Sistemi  
ABS Plastic Bottom Chassis and Central Brake System



5

Katlanabilir İtme Barı  
Foldable Push Bar



6

Sırt Kısımında Çift Amortisör  
Sırt Kısımında Çift Amortisör



**T-700  
HYDRAULIC  
EMERGENCY  
STRETCHER**



**FOOTREST  
POSITION**



**BACKREST  
POSITION**



**HEIGHT ADJUST.  
POSITION**



**TRENDELENBURG  
POSITION**



**R. TRENDELENBURG  
POSITION**

**T-700  
HYDRAULIC  
EMERGENCY  
STRETCHER**



## General Features

| Width<br>(Over all) | Lenght<br>(Over all) | Height<br>Min    | Height<br>Max    | Lying<br>Surface | Backrest<br>Angle | Legrest<br>Angle | Trendelenburg<br>Angle | R.trendelenburg<br>Angle | Safety<br>Working<br>Capacity |
|---------------------|----------------------|------------------|------------------|------------------|-------------------|------------------|------------------------|--------------------------|-------------------------------|
| <b>84<br/>cm</b>    | <b>208,5<br/>cm</b>  | <b>58<br/>cm</b> | <b>86<br/>cm</b> | 60x193 cm        | <b>90</b>         | <b>30</b>        | <b>10</b>              | <b>10</b>                | <b>230<br/>kgs</b>            |